



Unpack scheme:

1. Open the cardboard box.
2. Take out the can.
3. Release the tension ring lever and remove the top of the can.
4. Take out the upper polystyrene insert and the lead container.
5. Open the lead container and take out the vial holder.
6. Turn off the holder head and take out the vial with radiation material.

Package description:

Cardboard box dimensions: 220x220x230 mm

Tin can dimensions: 99x122 mm

Transport containers types: KT1-5 ÷ KT1-40

Weight: 1.5 ÷ 8.5 kg